


CONGJING MOXIBUSTION MODULE

CJ 360-1 ULTRA-SHORT-WAVE
THERAPEUTIC APPARATUS


Systemic treatment, dredge meridians, enhance immunity, bid farewell to the era of physical heat conduction, into the era of biological self-heating



OVERALL

In the Book of Rites, The Doctrine of the Mean, it goes: " In the world is the great foundation, and the way is the world. To neutralization, the status of heaven, all things yan."The evaluation and treatment of health and disease follow the path of neutralization. Yin Pingyang secret, qi and blood and, feeling calm, conform to heaven and earth for flat people,otherwise is not disease, disease, has been sick. And the disease is to have the imbalance of qi and blood, Yin and Yang, cold and heat, emotions and other conditioning to balance and state can be.


Physical heat: radiant heat and transfer



Physical Heat

- Radiant heat is a heating object that does not need to contact with the object that receives heat, and radiates heat to the receiving object in the form of waves.
- Conducheat means that the heating object contacts with the object that receives the heat and directly transmits the heat energy to the past. Physical heat has no energy form, is a direct transfer of heat.

Biological heat



Biological Heat


- Organisms under the RF field effect, make the body Electrolyte ions and charged particles for motion, and produce heat energy, is the biological energy into kinetic energy or potential energy, converted into heat, can the process, its advantage is from the inside out heat, make heat more uniform, not easy to burn, and makes the heat deeper, direct effect on the meridian qi and blood, is a safer and more high quality way of heat.

The action principle of the veiw-shot moxibustion chamber

- Microchip control technology and RF field effect technology are adopted to evenly distribute the RF energy in the body, so that the electrolyte ions and charged particles in the body can move, and then produce heat energy, namely the RF thermal effect. When the RF energy is low, the electrolyte ions and charged particles still move, producing some biophysical and biochemical reactions, but not the obvious thermal sensation, that is, the RF non-thermal effect.
- Radio frequency moxibustion therapy through the combination of radio frequency technology and traditional Chinese medicine treatment theory, effect on pain, hole or pathological meridians, promote meridian line the function of tissues and organs, cold clearing damp, through stasis, regulate qi, nourishing, adjust the balance of Yin and Yang, improve blood circulation, strengthen tissue nutrition, enhance immunity and phagocyte number, promote edema and inflammatory infiltration dissipate absorption, by reducing the nerve excitation, improve feeling threshold, achieve calm and alleviate the effect of muscle spasm.

The key technical specification

Service frequency	40.68MHz (Tolerance $\pm 1.5\%$)	Treatment timebased conditioning range	0 to 60 min, with error no greater than ± 1 min
Output power	800VA $\pm 20\%$. The 0~800VA is continuously adjustable	Display function	It has the real-time display function of incident power, reflection power and working time
Method of electrode output	Ring-shaped resonator (patented technology)	Sound cue function	With preheating and end of therapy
Resonator chamber body diameter	500mm	Indicator light prompt	Have various working status with the corresponding indicator light prompt
Resonant chamber body length	1300mm	Emergency shutdown function	The output power stops the output immediately after pressing the emergency stop button
Form size of the whole machine	Length: 2790mm width: 930mm height: 1245mm	Bed body in and out of the warehouse way	CJ 360-1 is electric / CJ 360-11 is manual



It can be adjusted according to the needs of different groups of people. Section output power and treatment time.


Personalized scheme

Treatment can act on several times other than the head. All areas, the treatment area is extensive.

Wide range of treatments



Product Characteristic



Ring speech vibration cavity
Crime technology utility model profit, use Reasonable matching, uniform distribution of the lifting field.

Two modes are optional
The treatment bed forward, back can be manually electric two operating modes are optional.



Scope of application

- Shoulder-hip soft tissue inflammation, pelvic tissue inflammation, sports injury and other local soft tissue inflammation, etc.
- Deficiency of qi and blood, imbalance of Yin and Yang, wind, cold and dampness, qi and blood stasis, phlegm dampness and obesity caused by meridian obstruction, liver deficiency, insomnia caused by heart and kidney failure, etc.



NMPA